
BODY FAT SCALES

MODEL - BS-04

INSTRUCTION MANUAL

BEDIENUNGSANLEITUNG
MODE D'EMPLOI
MANUALE DI ISTRUZIONI
MANUAL DE INSTRUCCIONES

Body Fat Scales - BS-04

PRODUCT SPECIFICATION


Product Dimensions: 300x300x20mm
11.8"x11.8"x0.8"

Power: 3xAAA Battery
Display: LCD

Max Weight: 180kg/400lbs/28st
Minimum Effective Weight: 6kg/13lbs

PRODUCT FEATURES

1. LCD display
2. Tempered glass surface
3. Battery Cover
4. Unit conversion: KG/LB/ST
5. High precision sensors
6. Metal sensor plates

SETTING UP THE PARAMETERS

The scales can store data on 10 separate users. To set up a new user or reset an existing user follow the steps below:

1. When the LCD display shows '0.0', press the 'SET' button for 3 seconds to enter the parameter input mode.
2. When the user code is flashing, press the up or down buttons to set the user number. Once selected press 'SET' to confirm your choice.
3. When the gender sign is flashing, press the up or down buttons to select the required gender. Once selected press 'SET' to confirm your choice.
4. Follow the same process for the height and age.

Remember your user code (1-10) to avoid setting the parameters again. If the scales are idle for more than 6 seconds during the parameter setting, the parameters already set will be saved and the LCD display will show "0.0" meaning that body fat measuring has now started.

WEIGHING

1. Place the scales on a solid and level surface.
2. Stepping onto the scales will automatically turn them on. Once a steady measurement is taken, the final weight measurement will flash on the screen.
3. Scales will automatically switch off when not in use.

MEASURING BODY FAT

For an accurate result, please ensure that shoes and socks are removed before stepping onto the unit. Then align your feet so that they are both touching the metal sensor plates.

1. Press the 'SET' button to switch the unit on.
2. Press the up or down button to select the user or reset the parameters.
3. When the correct user number is flashing, step onto the scales (ensuring both feet are touching the sensor plates).
4. Wait until the weight reading is locked and flashing. Remain on the scales while it calculates the other measurements.
5. Once the results have been calculated they will be displayed and cycled through. The results will repeat three times and saved for next time.

Body Fat %

Visceral Fat

Body Water %

Recommended Calorie Intake

Muscle Mass

Body Mass Index

Bone Mass

To display a user's last test results, follow the steps below:

1. Press 'SET' button to turn the scales on.
2. Select the correct user by pressing the up and down arrows.
3. Wait for 8 seconds until the user stops flashing or press firmly on the scales platform.
4. Press up and down arrows to cycle through last results.

AUTO IDENTIFY USER FUNCTION

After the user has set their parameters and weighed themselves, the scales will remember the user's weight and automatically revert to their user profile, if their weight is within +/- 2kg / 4.4lbs.

BODY MEASUREMENTS EXPLAINED

These bathroom scales use bioelectrical impedance to measure body fat. A very small and safe electrical current is passed through the lower half of the body. The electrical current flows slower through bone and fat than it does through muscle and water. The scales can measure this difference in speed and estimate body fat using a multi-step, mathematical formula.

Please note, these scales have been calibrated for those individuals who have a median body type and therefore the measurements recorded by individuals with unique body types e.g. significant upper body muscle, athletes or bodybuilders may not be accurate.

Body Fat %: This measurement refers to the subcutaneous fat (outer fat) closest to the skin. Men and women have a different distribution of body fat so it is normal for body fat % to vary based on gender.

AGE	MEN	WOMEN
20-39	8%-20%	22%-33%
40-59	11%-22%	24%-34%
60-79	13%-24%	25%-36%

Body Water %: Indicates the percentage of water in the body. It is important to stay well hydrated for overall health and wellbeing.

MEN	WOMEN
50%-65%	45%-60%

Muscle Mass: Indicates the total weight of muscle in your body. The normal ranges for muscle mass are as follows:

AGE	MEN	WOMEN
20-39	34%-40%	25%-30%
40-59	33%-39%	24%-30%
60-79	33%-38%	24%-30%

Bone Mass: Indicates the total weight of bone in your body. The normal ranges for bone mass are as follows:

MEN	WOMEN
3%-5%	2.5%-4%

Visceral Fat: Indicates the amount of visceral fat in your body. Visceral fat is found in the abdomen and surrounding vital organs.

Category	Visceral Fat
Normal	1-9
High	≥10

Recommended Calorie Intake: This calculates how many calories you should consume per day based on your data and body composition.

Body Mass Index: Body Mass Index (BMI) is a number calculated from a person's weight and height. Refer to the following chart for classifications:

Category	BMI
Underweight	< 18.5
Healthy weight	18.5-24.9
Overweight	25-29.9
Obese	30-34.9
Severely Obese	35-39.9
Morbidly Obese	> 40

NOTES

- Do not jump onto the scales as this may damage the measuring instruments.
- Always use on a hard surface. Measurements will be inaccurate on a soft surface such as carpet.
- Try and keep balanced and as steady as possible when standing on the scales.
- Always remove shoes and socks before using the scale.
- The body fat measurements are only optimised for people who are aged between 11-69.
- These scales are for home use and should not be used as a medical device.
- To get accurate results, it is recommended to weigh yourself at the same time of the day as measurements may vary after exercise, dieting or dehydration.

TO INSERT OR REPLACE BATTERY

1. Carefully remove the battery cover on the underside of the scales.
2. Insert 3 x AAA battery into the compartment as per the directions indicated.
3. After inserting the batteries, replace the battery cover.

UNIT CONVERSION

These scales measure weight in Kilograms (kg), Pounds (lb) and Stone (st). Press the button on the underside of the scales to switch between them.

ERROR CODES

- 'Lo' – Batteries are low and needs replacing.
- 'Err' – Scales are overloaded and have exceeded maximum weight. Step off the scales to avoid damaging the sensors.
- 'C' – A measuring error has occurred. Step off scales and back on again.

SAFETY

- Never overload the scales. This may damage the sensors.
- Do not leave heavy items on the scales when not in use.
- Keep away from water.
- To avoid slipping, keep the scales dry and never use with wet feet.
- Always keep the scales horizontal.
- To clean the scales, wipe with a soft damp cloth. Never use cleaning chemicals.
- The scales are is not recommended for the following:
 - o Pregnant women
 - o Those fitted with a pace maker or similar implanted medical device
 - o Those who have symptoms of edema
 - o Those of dialysis treatment
 - o Athletic or muscular body type (body builders)

DISPOSAL

This marking indicates that this product should not be disposed of with other household waste throughout the EU. This is to prevent possible harm to the environment and human wellbeing. Contact your local authority for information regarding the collection systems available. If possible, recycle it responsibly to promote the sustainable reuse of resources.

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© Copyright 2018